


# INSTRUCTION FOR USE


## Wheelchair table

Suitable models: H011 / H011B



- 1.TAKE THE TABLE OUT OF THE BOX
- 2.SCREW BOTH LEGS INTO THE TABLE
- 3.INSERT TABLE WITH LEGS WITH THE CLAMPS
- 4.ATTACH BOTH CLAMPS TO THE ARMRESTS OF THE WHEELCHAIR





Timago International Group  
Spółka z o.o. i Spółka – Spółka komandytowa  
ul. Karpacka 24/12  
43-316 Bielsko-Biała, Polska

T.: +48 33 499 50 00  
F.: +48 33 499 50 11  
E.: [info@timago.com](mailto:info@timago.com)



11/2020(I)